

The HF Mobile X-ray Generator consists of an inverter and oil-filled tank which houses the high voltage transformer and tube. It has been designed and developed for OEM application in Radiography for veterinary application. The high radiation output allows for lower kV per exposure, shorter exposure and superior imaging.

#### Features:

- ❖ Compact & light weight
- ❖ Reduction in dose
- ❖ Superior Imaging
- ❖ Tablet / Keyboard /touch-Screen Control
- ❖ Easy Serviceability
- ❖ OEM Customization Available



| TECHNICAL SPECIFICATIONS                  |   |  |   |
|---|---|--|---|
| PARAMETER                                 | SPECIFICATION   |  |   |
| <b>Input Voltage</b>                      | 230V $\pm$ 10% AC, 50Hz, single phase   |  |   |
| <b>Maximum Output Voltage &amp; Power</b> | <b>XRV5(5kW)</b><br>50mA at 100kV<br>100mA at 50kV  | <b>XRV10(10kW)</b><br>80mA at 120kV<br>250mA at 40kV | <b>XRV15(15kW)</b><br>125mA at 120kV<br>300mA at 50kV |
| <b>HF Generator Frequency</b>             | 80kHz   |  |   |
| <b>Ripple at the Maximum Power</b>        | <2%   |  |   |
| <b>Rise Time to Maximum Power</b>         | <3ms  |  |   |
| <b>Protections</b>                        | Over Voltage, Over Current, Over Temperature, Filament Fault, Earth Fault, Rotor Fault, Inverter Fault. |  |   |
| <b>Application</b>                        | Intensive care units and radiology department   |  |   |
| <b>kV Range</b>                           | 40kV to 100kV for 5kW<br>40kV to 120kV for 10kW & 15kW  |  |   |
| <b>mA Range</b>                           | 20mA to 300mA   |  |   |
| <b>Modes</b>                              | APR & Manual  |  |   |
| <b>mAs</b>                                | 0.1 to 50mAs  |  |   |
| <b>Exposure Time</b>                      | 4ms to 5sec   |  |   |
| <b>User Interface</b>                     | Keyboard with microcontroller-based control   |  |   |

For any queries or customization requests contact us at [info@ionics.co.in](mailto:info@ionics.co.in)

For product line information visit us at [www.ionics.co.in](http://www.ionics.co.in)